
## READ AND UNDERSTAND COMPLETE INSTRUCTIONS PRIOR TO USE

1. Remove product from packaging and inspect. If packaging is already open, or a component is missing or damaged, do not use.
2. Set up Bubble CPAP system or Nasal CPAP system.
3. Connect the humidification chamber and gas supply line to the outlet/inlet ports of the Pressure Relief Manifold.
4. Connect gas supply line to gas source.
5. If oxygen sensor is to be used, remove cap from port and insert sensor. Otherwise ensure cap is fully secured to oxygen sensor port.
6. Follow the manufacturer's instructions for use of the humidifier to ensure effective transfer of humid gas.

### Product Description

Pressure Relief Manifold with Gas Supply Tube and Oxygen Sensor Port is a single use device for use with patients weighing  $\leq 10$  kg. The pressure relief manifold is designed to prevent breathing circuit pressure from exceeding the preset pressure limit. Each manifold includes one 12.5 cm H<sub>2</sub>O pressure relief valve, gas supply line, and capped oxygen sensor port. The 22 mm ID inlet is compatible with all humidifier chambers.

It is designed to pair with Bubble CPAP and Nasal CPAP systems and can be used in neonatal critical care units, delivery rooms, pediatric critical care units, and special procedure units.

### Intended Use

The pressure relief manifold is designed to prevent breathing circuit pressure from exceeding the preset pressure limit. It is designed for use with patients weighing less than 10 kg.

### Storage

Store the device in a cool dry place out of direct sunlight before use.

### Disposal

Dispose of in accordance with local regulations.

Medical Technologies®  
www.bandb-medical.com

## Pressure Relief Manifold™

### 12.5 cm H<sub>2</sub>O Pressure Relief Manifold

Contents: 1 unit

**REF** 25330

**LOT** MO11030001

-20°C  
-4°F

Not Made With  
Natural Rubber Latex

**! READ PACKAGE  
INSERT FOR  
INSTRUCTIONS**

**USE BY:** 04/20/2027

**2** — Disposable —  
Single Patient Use Only  
**DO NOT REUSE**

**CAUTION:** Federal law  
restricts this device to sale by  
or on order of a physician.

**NON  
STERILE** — Packaged Clean  
— Not Sterile —

**Distributed By:**  
B&B Medical Technologies  
1954 Kellogg Ave. Ste 100  
Carlsbad, CA 92008 U.S.A  
Phone: +1.760.929.9972  
Fax: +1.760.929.9953  
www.bandb-medical.com

**EC REP** **EU Representative:**  
Mdi Europa GmbH  
Langenhagener Str. 71  
Germany  
Phone: +49 511 39089530  
Fax: +49 511 39089539

Made in Taiwan  
PN S533U Rev B 03/22/22  
© 2022 All Rights Reserved

### Caution

- This product should be used by trained personnel.
- Do not store in direct sunlight, extreme temperature or in high humidity. Store in dry conditions.
- Do not use if package is open or damaged.
- If oxygen sensor is not installed, ensure oxygen sensor port is sealed.

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