



Bite Proof Bite Block™

**Provides unobstructed access to oral cavity;
perfect for short term intubations, procedures**

Bite Proof Bite Block prevents the patient from biting the endotracheal tube and cuff inflating tube. The latex free Bite Proof Bite Block fits onto 5 – 9 mm endotracheal tubes and typically is used for short term intubations and procedures. The unique design allows for unobstructed access to the oral cavity. The Bite Proof Bite Block is placed on the endotracheal tube and the hypoallergenic adhesive-backed wing is placed on the cheek to secure the bite block placement. It can also be used for protection of closed suction catheters, endoscopes, bronchoscopes and oral gastric tubes.

When utilized with B&B's StabilTube™, LockTite™ and ET Tapes™, or with other products and adhesives used to secure and support the artificial airway, the Bite Proof Bite Block provides a fast and secure system for protecting the artificial airway.



SAFE

Latex free

Hypoallergenic materials for fragile skin

Provides convenient access for oral care, suctioning and procedures

Well tolerated by patients

COST EFFECTIVE

Choice for short term intubations and procedures

Packaged with securing device

CONVENIENT

Easy and quick application for time sensitive conditions

Adhesive tab ensures quick and easy application

VERSATILE

Fits all endotracheal tubes sizes 5 – 9 mm

Designed for use in critical care, surgery, emergency, transport and alternate care environments

Can be used in combination with B&B endotracheal tube holders and other commercially available securing devices

Ordering Information

11040 Bite Proof Bite Block (100/case)

11045 Bite Proof Bite Block (50/box)

Complementary Products

11060 StabilTube (40/case)

11010 ET Tape II (100/case)

11020 ET Tape (100/case)

11090 LockTite (50/case)

B&B Products are available from

B&B Medical Technologies and finer specialty care distributors.

Visit www.BandB-Medical.com or contact us today at 1.800.242.8778.

11040C2

© 2009 B&B Medical Technologies. All rights reserved.

"A World of Products for Better Breathing" is a service mark of B&B Medical Technologies.

Toll-free: +1.800.242.8778

Tel: +1.760.929.9972

Fax: +1.760.929.9953

2734 Loker Avenue West, Suite M
Carlsbad, CA 92010 USA

