



LockTite™

Secures ET tubes without adhesive tapes

LockTite is a patented, adjustable, non-adhesive endotracheal tube holder designed for use with patients who cannot tolerate adhesive products. It is ideal for patients with facial burns, severe facial irritation, beards or allergies to adhesives. LockTite is latex free and made of hypoallergenic materials. The pre-cut, non-irritating foam endotracheal tube holder is secured around the patient's neck with a VELCRO® brand neck band. A Releasable Cable Tie ensures the endotracheal tube stability and maximizes patient comfort while minimizing the incidence of inadvertent extubations.

B&B recommends using Universal Bite Block™ or Bite Proof Bite Block™ with LockTite to protect teeth and prevent biting on the endotracheal tube. LockTite can be used in all critical care units, operating rooms, emergency departments, special procedure units and transport situations.

Ordering Information

11090	LockTite (50/case)
-------	--------------------

Complementary Products

11160	Universal Bite Block (Adult) (50/case)
-------	--

11170	Universal Bite Block (Pediatric) (50/case)
-------	--

11040	Bite Proof Bite Block (100/case)
-------	----------------------------------

B&B Products are available from

B&B Medical Technologies and finer specialty care distributors.

Visit www.BandB-Medical.com or contact us today at 1.800.242.8778.



SAFE

Latex free, non-adhesive ET holder

Minimizes skin breakdown

Hypoallergenic materials for fragile skin

COST EFFECTIVE

Prepackaged kit with needed tools for securing the endotracheal tube

Reduces application time with VELCRO® brand fastener

Complete package for immediate use includes foam retaining unit and Releasable Cable Tie

CONVENIENT

Easy and quick application for time sensitive conditions

Adjustable fastener for quick application and change of position

VERSATILE

Used for endotracheal tubes 4 mm and greater

Child through adult application

Oral or nasal application

Designed for use in critical care and special procedure units, burn units and emergency situations

Accommodates patients weighing from 40 to 300 pounds

11090C2

© 2009 B&B Medical Technologies. All rights reserved. VELCRO® brand is a Registered Trademark of VELCRO Industries B.V. "A World of Products for Better Breathing" is a service mark of B&B Medical Technologies.

Toll-free: +1.800.242.8778

Tel: +1.760.929.9972

Fax: +1.760.929.9953

2734 Loker Avenue West, Suite M
Carlsbad, CA 92010 USA

B&B
MEDICAL TECHNOLOGIES
A World of Products for Better Breathing™