

Hybrid Baby Tape™

Latex Free Infant Endotracheal Tube Holder,
Cannula and NG/OG Tube Security
Hybrid Baby Tape (Part Number 11700)



Babies present with difficult endotracheal tube securing issues because of their potential for copious secretions, delicate skin integrity and frequent changes in position. These often adversely affect adhesive securing ability. Clinicians use configurations and cutting patterns of many different types attempting to stabilize the baby's endotracheal tube without harming tender skin. B&B Medical Technologies' Hybrid Baby Tape has the solution.

Hybrid Baby Tape combines the best features of hydrocolloid and adhesive tape. Gentle hydrocolloid tape safely secures and contours to the baby's tender skin while the adhesive tape firmly attaches to the endotracheal tube. The combination of tapes provides the perfect solution for long term attachment of the Hybrid Baby Tape and also for NasoGastric Tubes, OroGastric Tubes and Nasal Cannulas.

Benefits

Convenient...

- Easy and Quick Applications
- Time Saving, Pre Cut Tapes
- Complete Kit

Safe...



- Skin Protective Qualities
- Allows Skin to Breathe
- Latex Free, Hypo-Allergenic
- Enhanced Durability for Tube Protection

Cost Effective...

- Two Tapes per Package
- Secure Hold
- Used to Secure ETT, Cannulas and OG/NG tubes

Versatile...

- Infant and Pediatric Application
- Contours to Skin for Secure Fit
- Secures Additional Small Tubes and Cannulas
- Use in All Clinical Conditions



Ordering Information

11700 - Hybrid Baby Tape™ (100/Case)
11750 - Hybrid Baby Tape™ (50/Box)
Recommended for 2.0 to 3.5 mm ETTs

B&B Products are available from
B&B Medical Technologies
and Finer Specialty Care Distributors
Contact Us Today 1.800.242.8778
or www.BandB-Medical.com



B&B Medical Technologies
3800 Watt Avenue, Suite 275
Sacramento, CA 95821

Tel: +1.916.331.5221
Fax: +1.916.331.0161
www.BandB-Medical.com

A World of Products for Better Breathing™

B&B
MEDICAL TECHNOLOGIES