

Sil.Flex™ Stoma Pad

Revolutionary cushion absorbs pressure at stoma sites; multiple sizes for smallest infant to largest adult

Stoma sites require routine care throughout the day and may become sensitive or compromised due to constant pressure or movement of the tracheostomy tube and flange against surrounding tissue. Sil.Flex is a new ergonomically designed product that acts as a cushion between the flange and the stoma site.

The contoured surface of the Sil.Flex provides a stable, comfortable interface between the flange and the patient's neck. Sil.Flex Stoma Pads are made of a soft, fillerless elastomeric silicone encased in polyurethane. Hypoallergenic and latex-free, each sterile pad is individually packaged, and is intended for single patient application. Sil.Flex is easy to remove for patient hygiene and routine cleaning. The Sil.Flex Stoma Pad is available in four sizes, accommodating small infants to large adults.

The TrachGuard™ or the TrachStay™ may be used in combination with the Sil.Flex Stoma Pad. Sil.Flex Stoma Pads enhance patient comfort with all brands and styles of tracheostomy tubes.



The Sil.Flex Stoma Pad can be used in all critical care units, operating rooms, emergency departments and special procedure units as well as alternate care facilities and home care. Using Sil.Flex Stoma Pad may assist in reducing irritation at the stoma site.

SAFE

Contoured surface provides pressure absorbing cushion that may reduce irritation at the stoma site

Hypo-allergenic, latex free, medical grade silicone pads

Size clearly identified on each pad

Flat surface conforms to tracheostomy tube flange

Center opening reduces direct tube contact

COST EFFECTIVE

Individually packaged, sterile and easy to apply

28 day application, easily cleaned at patient bedside

CONVENIENT

Available in four sizes, accommodating small infants through large adults





Minimizes clinician time for application with individually packaged Sil.Flex Stoma Pad

VERSATILE

Unique sizes and shapes to fit different shaped necks and tracheal tube sizes

Special design allows for placement without tube removal

Ordering Information

20000	Sil.Flex Stoma Pad – Size 0 (10/box)	
20010	Sil.Flex Stoma Pad – Size 0 (40/case)	
20001	Sil.Flex Stoma Pad – Size 1 (10/box)	
20011	Sil.Flex Stoma Pad – Size 1 (40/case)	
20002	Sil.Flex Stoma Pad – Size 2 (10/box)	
20012	Sil.Flex Stoma Pad – Size 2 (40/case)	
20003	Sil.Flex Stoma Pad – Size 3 (10/box)	
20013	Sil.Flex Stoma Pad – Size 3 (40/case)	

Complementary Products

11250	TrachGuard
11270	TrachStay

B&B Products are available from

B&B Medical Technologies and finer specialty care distributors.

Visit www.BandB-Medical.com or contact us today at 1.800.242.8778.

Toll-free: +1.800.242.8778

Tel: +1.760.929.9972

Fax: +1.760.929.9953

2734 Loker Avenue West, Suite M
Carlsbad, CA 92010 USA



20000C3

© 2010 B&B Medical Technologies. All rights reserved.

Sil.Flex Stoma Pad is a trademark of A Plus Medical.

"A World of Products for Better Breathing" is a service mark of B&B Medical Technologies.