

Sat Point™ Pedi



B&B Sat Point™ Pedi Oximeter is specifically designed for infant and pediatric patients. The Sat Point Pedi measures SpO2, Pulse Rate and provides a clinical tool needed to identify changes in patient condition. The compact, comfortable finger design with adjustable brightness and color display provides flexibility for use in clinics, hospital, homecare, physician office and EMS environments. Most Important, the Sat Point™ Pedi has the best fit for the smallest of fingers. Sat Point™ Pedi: Part Number 22020

Benefits

Convenient...

Color Display of SpO2, HR, Waveform
Pulse Strength Signal Display
View Data in Multiple Displays

Safe...

Finger Out Alarm, Automatic Power Off
Pulse Rate Audio
Indicator to Identify Best Finger Placement

Technical Specifications

Patient Range: Infant & Pediatric
SpO2: Range: 0 to 99%/Resolution: 1%
Accuracy: 80 to 99%: ± 2%
70 to 80%: ± 3%
0 to 69%: unspecified
Heart Rate/Pulse:
Range: 0 to 254 bpm/Resolution: 1 bpm
Accuracy: ± 2bpm or ± 2%
Audible: PR Tone modulation
Information: Signal Strength, Bar Graph
Alarm: Finger Out, Low Power Indicator
Display Parameters: SpO2, PR, SpOs waveform

Cost Effective...

Lower Power Consumption with Low Battery Indicator
AAA Battery Power for up to 30 hours of operation
Two Year Warranty

Versatile...

Infant and Pediatric Application
Six Display Modes/Adjustable Brightness

Environmental:

Operating Temperature: 5C to 40 C
Storage Temperature: -10C to 40 C
Operating Humidity: 15% RH to 80% RH
Storage Humidity: 10% RH to 80% Rh

Mechanical:

Size: 58 x 32 x 28 mm / Weight: 50 gm (w/batteries)
Power: Type: 2 AAA Alkaline Batteries
Operating Time: 30 hours nominal use
Display: LCD/6 display modes
Brightness: adjustable levels/10
Please refer to the instructions for use with the Sat Point Pedi for detailed instructions, warnings and cautions.

Ordering Information

22020 - Sat Point™ Pedi Oximeter (1/box)

**B&B Products are available from
B&B Medical Technologies
and Finer Specialty Care Distributors
Contact Us Today 1.800.242.8778
or www.BandB-Medical.com**



B&B Medical Technologies
2734 Loker Ave West-Suite M
Carlsbad, CA 92010-6628

Tel: +1.760.929.9972
Fax: +1.760.929.9953
www.BandB-Medical.com

A World of Products for Better BreathingSM

B&B
MEDICAL TECHNOLOGIES