

Sat Point™ Plus



B&B Sat Point™ Plus Oximeter is the advanced finger model for pediatric and adult patients. The Sat Point Plus measures SpO2, Pulse Rate and provides a clinical tool needed to identify changes in patient condition. The compact, comfortable finger design with adjustable brightness and multi-color display provides flexibility for use in clinics, hospital, homecare, physician office and EMS environments. Most Important, the Sat Point Plus™ can cover a broad range of patients.
Sat Point™ Plus: Part Number 22030

Benefits

Convenient...

MultiColor Display of SpO2, HR, Waveform
Pulse Strength Signal Display
View Patient Data in Multiple Displays

Safe...

Finger Out Alarm, Automatic Power Off
Pulse Rate Audio

Technical Specifications

Patient Range: Adult & Pediatric
SpO2: Range: 0 to 99%/Resolution: 1%
Accuracy: 80 to 99%: $\pm 2\%$
70 to 80%: $\pm 3\%$
0 to 69%: unspecified
Heart Rate/Pulse:
Range: 0 to 254 bpm/Resolution: 1 bpm
Accuracy: ± 2 bpm or $\pm 2\%$
Audible: PR Tone modulation
Information: Signal Strength, Bar Graph
Alarm: Finger Out, Low Power Indicator
Display Parameters: SpO2, PR, SpOs waveform

Cost Effective...

Lower Power Consumption with Low Battery Indicator
AAA Battery Power for up to 30 Hours of Operation
Two Year Warranty

Versatile...

Adult and Pediatric Application
Six Display Modes/Adjustable Brightness

Environmental:

Operating Temperature: 5C to 40 C
Storage Temperature: -10C to 40 C
Operating Humidity: 15% RH to 80% RH
Storage Humidity: 10% RH to 80% Rh

Mechanical:

Size: 58 x 32 x 28 mm / Weight: 50 gm (w/batteries)
Power: Type: 2 AAA Alkaline Batteries
Operating Time: 30 hours nominal use
Display: LCD/6 display modes
Brightness: adjustable levels/10
Please refer to the instructions for use with the Sat Point Plus for detailed instructions, warnings and cautions.

Ordering Information

22030 - Sat Point™ Plus Oximeter (1/box)

**B&B Products are available from
B&B Medical Technologies
and Finer Specialty Care Distributors
Contact Us Today 1.800.242.8778
or www.BandB-Medical.com**



B&B Medical Technologies Tel: +1.760.929.9972
2734 Loker Avenue West-Suite M Fax: +1.760.929.9953
Carlsbad, CA 92010-6628 www.BandB-Medical.com

A World of Products for Better BreathingSM

B&B
MEDICAL TECHNOLOGIES