



Precision Resistor Kit™

Ideal teaching tool demonstrates, simulates changes in resistance

Precision Resistor Kit provides the healthcare practitioner with calibrated adapters to demonstrate changes in airway resistance. The Precision Resistor Kit contains three resistors: Rp5, Rp20 and Rp50 and simulates nominal levels of resistance. Used in tandem with The Test Lung™, The Precision Resistor Kit is the ideal tool for teaching and demonstration of changes in resistance.

Ordering Information

20118	Precision Resistor Kit (1/box)
-------	--------------------------------

Complementary Products

20110	The Test Lung 1.0 L (Adult) (1/box)
-------	-------------------------------------

20120	The Test Lung 0.5 L (Pediatric) (1/box)
-------	---

SAFE

Factory calibrated

Easily cleaned and sterilized between uses

COST EFFECTIVE

Complete package of three resistors: Rp5, Rp20, Rp50

CONVENIENT

Easy to use

Quick application

VERSATILE

Designed for use with demonstration and testing of ventilators in respiratory care, clinical teaching areas and biomedical departments



The Test Lung 1.0 L (Adult)



The Test Lung 0.5 L (Pediatric)

B&B Products are available from

B&B Medical Technologies and finer specialty care distributors.

Visit www.BandB-Medical.com or contact us today at 1.800.242.8778.

20118C1

© 2009 B&B Medical Technologies. All rights reserved.
 "A World of Products for Better Breathing" is a service mark of B&B Medical Technologies.

Toll-free: +1.800.242.8778
 Tel: +1.760.929.9972
 Fax: +1.760.929.9953
 2734 Loker Avenue West, Suite M
 Carlsbad, CA 92010 USA

