

Single Pole Clamp



Dual Pole Clamp



# Venti.Plus™ Single Pole Clamp Venti.Plus Dual Pole Clamp

Unique clamps for positive attachment to standard equipment poles

The Venti.Plus Single Pole Clamp and Venti.Plus Dual Pole Clamp are sturdy pole clamps designed to work in tandem with the Babi.Plus™ Bubble PAP Valve 0–10 cm H<sub>2</sub>O. The Venti.Plus Single Pole Clamp is intended to support up to 3 kg in weight, just right for the Babi.Plus Bubble PAP Valve. The Venti.Plus Dual Pole Clamp is well-built and intended to support up to 10 kg of weight, which conveniently allows for mounting of the Babi.Plus Bubble PAP Valve and a humidifier or other devices.

The Venti.Plus pole clamps are ISO compliant for mounting on poles 2.0 to 3.2 cm in diameter.

The Venti.Plus pole clamps are a complementary product of the Babi.Plus Bubble PAP Valve product line, and are available separately.



## SAFE

Specifically designed to hold Babi.Plus Bubble PAP Valve

The dual pole clamp can be used for multiple devices weighing ≤ 10 kg

## COST EFFECTIVE

Sturdy clamp that can be used with multiple devices

Unique dual pole clamp minimizes expense of two pole clamps

## CONVENIENT

Easy and quick application with thumb screw mechanism to standard equipment pole

Easily cleaned between uses

## VERSATILE

Designed for use with Babi.Plus Bubble PAP Valve, humidifier or other devices

## Ordering Information

20271	Venti.Plus Single Pole Clamp (1 each)
20272	Venti.Plus Dual Pole Clamp (1 each)

## Complementary Product

20205	Babi.Plus Bubble PAP Valve 0–10 cm H <sub>2</sub> O (12/box)
20210	Babi.Plus Bubble PAP Valve 0–10 cm H <sub>2</sub> O (24/case)

## B&B Products are available from

B&B Medical Technologies and finer specialty care distributors.  
Visit [www.BandB-Medical.com](http://www.BandB-Medical.com) or contact us today at 1.800.242.8778.

20271C2

© 2011 B&B Medical Technologies. All rights reserved.  
Babi.Plus and Venti.Plus are trademarks of A Plus Medical.  
"A World of Products for Better Breathing" is a service mark of B&B Medical Technologies.

Toll-free: +1.800.242.8778  
Tel: +1.760.929.9972  
Fax: +1.760.929.9953  
2734 Loker Avenue West, Suite M  
Carlsbad, CA 92010 USA

