



Accutach™ Pneumotach Flow Sensor

Accurate and reliable flow rate monitoring for pediatric and adult patients

The Accutach sensor is a latex free, hypo-allergenic, flow sensing device intended for single patient use. When placed in the line of the ventilator circuit the sensor accurately measures flow, volume, and pressures delivered to the patient. This disposable device provides differential pressure signals through the sensor ranges, precisely measuring ventilator settings and assuring the integrity of the patient circuit.

The Accutach operates at the distal end of the ventilator circuit, located between the patient Y and the endotracheal tube. The small dead space of 11 cc and the Accutach's light weight provide ideal patient monitoring. The Accutach comes as a complete package with a Flow Sensor, a pair of 72 inch color-coded ventilator panel lines, a 15 mm ID – 22 mm OD adaptor, and a 22 mm tubing clip.

Each device is packaged clean and must be calibrated prior to use. The Accutach was designed to work with the H7200 VR-1 monitor and the Galileo®, Amadeus®, Veolar®, Raphael®, G5® and C-2® ventilators.

Ordering Information

25210	Accutach Flow Sensor (10/box)
25250	Accutach Flow Sensor (50/case)

Technical Parameters

Flow	0/min – 180/min
Dead space	11 cc
Resistance	< 1.6 mbar/1/s
Weight	11 g without tubing, 60 g total
Connectors	ISO 15 male – 15 male
Accuracy	+/- 25 ml or +/-10%, uncalibrated

B&B Products are available from

B&B Medical Technologies and finer specialty care distributors.
Visit www.BandB-Medical.com or contact us today at 1.800.242.8778.

SAFE

Latex free, hypo-allergenic materials

Accurate and reliable flow rate measurement

COST EFFECTIVE

Complete package with Flow Sensor, adaptor, ventilator panel lines, and tubing clip

An inexpensive method of obtaining ventilator-measured values

CONVENIENT

Color-coded ventilator panel lines for ease of application

Tubing Clip prevents inadvertent kinking

VERSATILE

Can be used for pediatric and adult patients

Designed to work with the H7200 VR-1 monitor and the Galileo, Amadeus, Veolar, Raphael, G-5, and C-2 ventilators

25210C2

© 2013 B&B Medical Technologies. All rights reserved.
Accutach is a trademark of Carlsbad International Export Inc.
Galileo, Amadeus, Veolar, Raphael, G5 and C-2 are registered trademarks of Hamilton Medical.



Toll-free: +1.800.242.8778
Tel: +1.760.929.9972
Fax: +1.760.929.9953
1954 Kellogg Avenue, Suite 100
Carlsbad, CA 92008 USA