



## ET Tape™

**Minimizes inadvertent extubations;  
ensures tube stability**

ET Tape is a pre-cut, non-irritating foam adhesive endotracheal tube holder designed to ensure stability of the endotracheal tube and minimize inadvertent extubations. ET Tape's latex free hypoallergenic materials have superior holding power and protect the upper lip and skin across the cheek area. A protective skin prep pad is provided to form a topical barrier between skin and adhesive for added protection against skin irritation and breakdown. A releasable cable tie secures placement of the endotracheal tube and prevents slippage. ET Tape provides tube stability for short-term applications and allows for consistency of taping and securement. ET Tape comes as a prepackaged kit that includes tape, skin prep pad and releasable cable tie.

B&B recommends using Universal Bite Block™ or Bite Proof Bite Block™ with ET Tape to protect teeth and prevent biting on the endotracheal tube. ET Tape can be used in all critical care units, operating rooms, emergency departments, special procedure units and transport situations.

### Ordering Information

---

11020	ET Tape (100/case)
-------	--------------------

---

### Complementary Products

---

11160	Universal Bite Block (Adult) (50/case)
-------	--

---

11170	Universal Bite Block (Pediatric) (50/case)
-------	--

---

11040	Bite Proof Bite Block (100/case)
-------	----------------------------------

---

### B&B Products are available from

B&B Medical Technologies and finer specialty care distributors.  
Visit [www.bandb-medical.com](http://www.bandb-medical.com) or contact us today at 1.800.242.8778.

#### SAFE

Latex free, hypoallergenic

Skin prep pad provides topical barrier to minimize skin irritation and breakdown

Enhanced durability for endotracheal tube protection and position

Non absorbent material reduces tube slippage

AHA guidelines recommend use of a commercial pre-cut device for ET tube security

#### COST EFFECTIVE

Prepackaged kit for immediate use includes pre-cut tape, skin prep and releasable cable tie

Reduces cost of repeated taping procedures

#### CONVENIENT

Pre-cut tape minimizes application time

Easy and quick application for time sensitive conditions

#### VERSATILE


Used for endotracheal tubes 4 mm and greater

Child through adult application

Oral or nasal application

Designed for use in critical care and special procedure units, operating rooms, emergency and transport situations

Accommodates patients weighing 40 to 300 pounds



11020C4

© 2014 B&B Medical Technologies. All rights reserved.

**BB**  
BREATHING INNOVATION



Medical Technologies  
[www.bandb-medical.com](http://www.bandb-medical.com)

Toll-free: +1.800.242.8778

Tel: +1.760.929.9972

Fax: +1.760.929.9953

1954 Kellogg Avenue, Suite 100

Carlsbad, CA 92008 USA