



# TrachGuard™

## Complete kit prevents inadvertent disconnects; stabilizes ventilator circuit and tracheostomy tube

TrachGuard is a patented combination long term tracheostomy tube holder and anti-disconnect device that dramatically reduces the potential for inadvertent ventilator circuit disconnects and exposure to cross contamination. Made to accommodate any tracheostomy tube, the anti-disconnect device stabilizes the patient connection to the ventilator circuit or in-line suction device.

TrachGuard is a complete kit that secures the tracheostomy tube to the patient connector of the ventilator circuit and a comfortable soft collar with hook and loop closure to secure the tracheostomy tube in position. The TrachGuard collar has superior holding power and is made of latex free, hypoallergenic soft foam.

The collar fits neck sizes 7 to 18 inches. TrachGuard is an essential tool for busy ICUs, long term acute care, home and alternate care facilities.



### SAFE

- Latex free, hypoallergenic materials
- Enhanced durability for securing the ventilator circuit connection
- Minimizes potential for decannulation
- Reduces potential exposure to contaminated tubing moisture

### COST EFFECTIVE

- Prepackaged kit with collar and anti-disconnect device
- Designed for use in critical care, PICU, long term acute care, alternate site and home care

### CONVENIENT

- Prepackaged kit minimizes time for setup and connection
- Easy and quick application to all tracheostomy tubes
- Comfortable foam collar fits pediatric through adult patients

### VERSATILE

- Adult and pediatric applications
- Collar fits neck sizes 7 to 18 inches
- Can be used with delivery of aerosol and supplemental oxygen therapy

### Ordering Information

11250	TrachGuard (80/case)
11255	TrachGuard (20/box)

### B&B Products are available from

B&B Medical Technologies and finer specialty care distributors. Visit [www.BandB-Medical.com](http://www.BandB-Medical.com) or contact us today at 1.800.242.8778.

