



# LockTite™

## Secures ET tubes without adhesive tapes

LockTite is a patented, adjustable, non-adhesive endotracheal tube holder designed for use with patients who cannot tolerate adhesive products. It is ideal for patients with facial burns, severe facial irritation, beards or allergies to adhesives. LockTite is latex free and made of hypoallergenic materials. The pre-cut, non-irritating foam endotracheal tube holder is secured around the patient's neck with a VELCRO® brand neck band. A Releasable Cable Tie ensures the endotracheal tube stability and maximizes patient comfort while minimizing the incidence of inadvertent extubations.

B&B recommends using Universal Bite Block™ or Bite Proof Bite Block™ with LockTite to protect teeth and prevent biting on the endotracheal tube. LockTite can be used in all critical care units, operating rooms, emergency departments, special procedure units and transport situations.

### SAFE

Latex free, non-adhesive ET holder  
Minimizes skin breakdown  
Hypoallergenic materials for fragile skin

### COST EFFECTIVE

Prepackaged kit with needed tools for securing the endotracheal tube  
Reduces application time with VELCRO® brand fastener  
Complete package for immediate use includes foam retaining unit and Releasable Cable Tie

### CONVENIENT

Easy and quick application for time sensitive conditions  
Adjustable fastener for quick application and change of position

### VERSATILE

Used for endotracheal tubes 4 mm and greater  
Child through adult application  
Oral or nasal application  
Designed for use in critical care and special procedure units, burn units and emergency situations  
Accommodates patients weighing from 40 to 300 pounds

### Ordering Information

11090	LockTite (50/case)
-------	--------------------

### Complementary Products

11160	Universal Bite Block (Adult) (50/case)
11170	Universal Bite Block (Pediatric) (50/case)
11040	Bite Proof Bite Block (100/case)

### B&B Products are available from

B&B Medical Technologies and finer specialty care distributors.  
Visit [www.bandb-medical.com](http://www.bandb-medical.com) or contact us today at 1.800.242.8778.



11090C3

© 2016 B&B Medical Technologies. All rights reserved. VELCRO® brand is a Registered Trademark of VELCRO Industries B.V.

**B|B**  
BREATHING INNOVATION



Medical Technologies  
[www.bandb-medical.com](http://www.bandb-medical.com)

Toll-free: +1.800.242.8778  
Tel: +1.760.929.9972  
Fax: +1.760.929.9953  
1954 Kellogg Avenue, Suite 100  
Carlsbad, CA 92008 USA