



HOPE™ Nebulizer

Delivers outstanding nebulizer performance with patented Heliox port

The HOPE Nebulizer has proven to be a cost effective solution for delivery of continuous medication for the treatment of moderate to severe asthma and compromised airways. The HOPE is a high output nebulizer used for delivery of aerosolized medications and diagnostic formulations. It provides additional hydration for patients to loosen secretions when extended therapy is required.

The HOPE is a refillable, closed dilution, high output, large volume continuous medication nebulizer. The flow design and 200 ml reservoir allow for up to 8 hours of nebulized medication delivery at 25 ml/hr. The HOPE delivers 1 to 8 hours of continuous medication with 2.5 MMAD. The output at 25 ml/hr delivers >60% particles in the respirable range at 10 lpm flow. With the addition of supplemental Heliox, the HOPE provides a 50 ml/hr output with >75% particles in the respirable range with 15 lpm supplemental 80/20 Heliox delivered through the auxiliary port.

The HOPE is the gold standard for delivery of continuous nebulization and Heliox therapy. Heliox has been reported to build a therapeutic bridge in severe asthma and compromised airways until conventional medications have had time to take effect. B&B's patented supplemental port for Heliox delivery greatly enhances the performance of the HOPE Nebulizer. It is no wonder that the competition has reconfigured their products to emulate the HOPE Nebulizer. B&B first introduced the patented supplemental port for Heliox delivery in a closed dilution nebulizer, greatly enhancing the performance of the HOPE Nebulizer. HOPE is the industry leader in simplicity and outstanding nebulizer performance in a Large Volume Continuous Medication Nebulizer.

Ordering Information

11310	HOPE Nebulizer (25/case)
11320	HOPE Nebulizer Adult Kit (10/case) Nebulizer, 6' Aerosol Tubing and Adult Mask
11330	HOPE Nebulizer Pediatric Kit (10/case) Nebulizer, 6' Aerosol Tubing and Pediatric Mask



HOPE Nebulizer Kit

B&B Products are available from

B&B Medical Technologies and finer specialty care distributors.
Visit www.bandb-medical.com or contact us today at 1.800.242.8778.

SAFE

Rescue therapy for asthmatics and obstructive lung disease patients

Latex-free solution

Easy-to-understand dosing chart

Allows clinician to add auxiliary gas flow without recalculation of medication

Hypoallergenic materials

COST EFFECTIVE

Prepackaged kit with needed accessories for delivery of continuous medication therapy

Dosing chart provided with each HOPE Nebulizer

Delivers medication therapy and humidification for up to 8 hours

CONVENIENT

Quick connection to flowmeter

Easy and quick application for time-sensitive conditions

Twist off reservoir for easy medication changes

Standard auxiliary port for delivery of therapeutic gases

VERSATILE

Designed for use in the critical care units, emergency department, extended care clinics, outpatient clinics and in hospital transport

Accommodates auxiliary port for delivery of supplemental heliox gas for therapy

11310C7

© 2021 B&B Medical Technologies. All rights reserved. The HOPE Nebulizer is a trademark of B&B Medical Technologies.

B|B
BREATHING INNOVATION



Medical Technologies®
www.bandb-medical.com

Toll-free: +1.800.242.8778

Tel: +1.760.929.9972

Fax: +1.760.929.9953

1954 Kellogg Avenue, Suite 100

Carlsbad, CA 92008 USA