



# StabiTube™

## Long-term protection eliminates repetitive taping; neck band adjusts for changing conditions

StabiTube is a patented endotracheal tube holder specifically designed to eliminate repetitive taping and is intended for use as an immediate endotracheal tube holder and stabilization device. Ease of use and quick application make StabiTube a favorite in time sensitive situations. StabiTube is a pre-cut, latex free, hypoallergenic foam adhesive ET holder that protects the upper lip and skin across the cheek area. Its adjustable, flexible VELCRO® brand neck band accommodates changes in patient's condition, positioning and neck size without needing to retape. StabiTube ensures ET tube stability and minimizes the incidence of inadvertent extubations. A Releasable Cable Tie secures the endotracheal tube and prevents slippage while also allowing for tube depth readjustment. A protective skin wipe is provided to form a topical barrier between skin and adhesive for added protection against skin breakdown and irritation. StabiTube comes as a complete kit that includes pre-cut tape, skin prep and Releasable Cable Tie.



B&B's Universal Bite Block™ or the Bite Proof Bite Block™ are recommended for use with StabiTube to protect the patient's teeth and prevent the patient from biting on the endotracheal tube. The StabiTube can be used in all critical care units, operating rooms, emergency departments, special procedure units and for transport situations.

### SAFE

Latex-free solution

Minimize skin breakdown

Hypoallergenic materials for fragile skin

Enhanced durability for endotracheal tube protection and position

### COST EFFECTIVE

Long-term solution for securing the endotracheal tube reduces the cost of repeated taping procedures

Complete package for immediate use includes tape, skin prep and Releasable Cable Tie

### CONVENIENT

Pre-cut tape minimizes clinician time for preparation and allows quick application for time sensitive conditions

Non absorbent material to reduce tube slippage

Easy adjustment with flexible VELCRO® brand neckband for change in patient condition

### VERSATILE

Used for endotracheal tubes 4 mm and greater

Child through adult application

Oral or nasal application

Designed for use in the all critical care units and transport situations

Accommodates patients from 40 to 300 pounds, neck sizes up to 26 inches (C-collar application)

### Ordering Information

11060	StabiTube (40/case)
-------	---------------------

11065	StabiTube (20/box)
-------	--------------------

### Complimentary Products

11160	Universal Bite Block (Adult) (50/case)
-------	--

11170	Universal Bite Block (Pediatric) (50/case)
-------	--

11040	Bite Proof Bite Block (100/case)
-------	----------------------------------

### B&B Products are available from

B&B Medical Technologies and finer specialty care distributors.

Visit [www.bandb-medical.com](http://www.bandb-medical.com) or contact us today at 1.800.242.8778.

11060C3

© 2016 B&B Medical Technologies. All rights reserved. VELCRO® brand is a Registered Trademark of VELCRO Industries B.V.

**B|B**  
BREATHING INNOVATION



Medical Technologies  
[www.bandb-medical.com](http://www.bandb-medical.com)

Toll-free: +1.800.242.8778

Tel: +1.760.929.9972

Fax: +1.760.929.9953

1954 Kellogg Avenue, Suite 100

Carlsbad, CA 92008 USA