

B&B Single Pole Clamp



B&B Single Pole Clamp

Unique clamp for positive attachment to standard equipment poles

The B&B Single Pole Clamp is a sturdy pole clamp that works in tandem with the B&B Bubbler Water Seal CPAP Valve 0–10 cm H₂O. The B&B Single Pole Clamp is intended to support up to 3.6 kg in weight, just right for the B&B Bubbler. The B&B Single Pole Clamp is ISO compliant for mounting on poles 2.0 to 3.2 cm in diameter.

The B&B Single Pole Clamp is a complementary product of the B&B Bubbler and is available separately.

SAFE

Compatible with the B&B Bubbler Water Seal CPAP Valve

COST EFFECTIVE

Sturdy clamp that can be used with multiple devices

CONVENIENT

Easy and quick application with knob tightening mechanism to standard equipment pole

Easily cleaned between uses

VERSATILE

Designed for use with the B&B Bubbler or other devices



B&B Bubbler

Ordering Information

| | |
|-------|--------------------------------|
| 25371 | B&B Single Pole Clamp (1 each) |
|-------|--------------------------------|

Complementary Product

| | |
|-------|---|
| 25312 | B&B Bubbler Water Seal CPAP Valve 0–10 cm H ₂ O (12/box) |
|-------|---|

| | |
|-------|---|
| 25324 | B&B Bubbler Water Seal CPAP Valve 0–10 cm H ₂ O (24/box) |
|-------|---|

B&B Products are available from

B&B Medical Technologies and finer specialty care distributors.

Visit www.bandb-medical.com or contact us today at 1.800.242.8778.