



Baby Tape Plus™

**Hydrocolloid Infant endotracheal tube holder
safely secures, contours to tender skin.**

Baby Tape Plus™ is a pre-cut endotracheal tube holder made from a non-irritating, latex free hydrocolloid adhesive that is skin friendly for infants. Hydrocolloid tape safely secures and contours to the baby's tender skin while the adhesive tape firmly attaches to the endotracheal tube. Baby Tape Plus™ is designed to be quick and simple to use in a wide variety of situations. The circular tube wrap makes for a cleaner, easier application and the larger wing tabs cover a bigger surface area for greater security, they can be trimmed to fit a wide range of infants. The middle tape strips can also be used to secure nasogastric /orogastric tubes and cannulas.

SAFE

Hydrocolloid adhesive is gentle on babies' delicate skin

Minimizes skin reaction

Hypoallergenic, non-irritating, latex free materials

COST EFFECTIVE

Prepackaged tapes reduce time to secure endotracheal tubes

Individual packages contain two tapes each

Stays in place during routine patient movement and transport

CONVENIENT

Easy and quick application for time sensitive conditions

Circular tube wrap makes for a cleaner, easier application

Larger wing tabs cover a bigger surface area for greater security

VERSATILE

Recommended for use with patients ≤ 10kg

Wing Tabs can be trimmed to fit a wide range of infants

One product can secure endotracheal tubes, nasogastric /orogastric tubes and cannulas

Ordering Information

25390	Baby Tape Plus (100/case - Contains 50 packs of 2 each)
-------	---

25395	Baby Tape Plus (50/box - Contains 25 packs of 2 each)
-------	---

Complementary Products

11100	NeO2-Safe (100/case)
-------	----------------------

B&B Products are available from

B&B Medical Technologies and finer specialty care distributors.

Visit www.bandb-medical.com or contact us today at 1.800.242.8778.

